

L.E.D. CURING LIGHTS



FX 35-P



FX 35



FX 45

FX 35 / 35-P Power Pen

- 8mm Tip with Anti-glare Shield
- Powerful L.E.D. Monitor
- Start/Stop Trigger
- Sterilizable Curing Tip Protects L.E.D. Module
- Pre-select Timer- 10, 20 or 30 sec.

The **FX 35-P** includes the same features as the FX 35 plus a built-in fan providing constant cool down and allowing for extended use.

FX 45 Cordless L.E.D. Curing Light

- Powerful L.E.D. Monitor
- Fast Cure Cycle as well as Ramp, Constant & Pulse Curing Cycles
- Built-in Radiometer
- Digital Display of Curing Mode & Timer
- Long Lasting & Replaceable Li-ion Battery
- Up to 120 min. Run Time per Charge
- Sterilizable Curing Tip Protects L.E.D. Module
- Battery Charging Station with Light Intensity Meter